



A DYNAMIC KNEE EXAM



Evidence-based care¹



Patient satisfaction
and positive outcomes²



Value to patient
and practice²

A Technological Edge in
Decision-Making to better
**CORRECT KNEE
ALIGNMENT AND
RESTORE 3D
KNEE FUNCTION**






THE FIRST

Technology to Accurately Measure 3D Dynamic Alignment that Delivers Objective Data


(FDA Cleared, Health Canada and CE Marked Class II Medical Device)*



ACCURACY Now at your Fingertips!

KneeKG, proven¹ high accuracy to assess 3D function

KneeKG lower-limb alignment correlates with alignments from both robotic assisted system and long leg standing X-ray⁵



Increasing EFFICIENCY AND EFFICACY OF YOUR TOTAL JOINT PROGRAM for Surgical and Non-Surgical Candidates

3D Dynamic Malalignment is Linked to:

- Poor patient reported symptoms⁶
- Poor outcomes (surgical and non-surgical interventions)⁷
- Bone marrow lesions⁸
- Patient dissatisfaction⁷

- 5Y Contribution Margin: **\$1M-\$2M**
- Y1 #Uses: **650**; Yearly growth: **10%**
- Reimb: **\$220-\$400** > Range based upon payer reimbursement rates reported by customers

Back to performance

"KneeKG helped me better understand my knee pain and generate a personalized treatment plan, which allowed me to be the same athlete I was before my injury."

- Alessandro Riggi
Professional Soccer Player

Knee osteoarthritis

"I changed my standard exercises for the specific plan recommended, my knee pain went from 10/10 to 0 and I am now back to my activities"

- Bruny Surin
Olympic Gold Medalist, 4x100m



Painful knee post surgery

"KneeKG started me down a whole new path!"

It was so good to finally learn why I was experiencing pain and have an action plan. If only I had access to KneeKG long ago, I believe I could have avoided many painful months and been so much further along in my recovery."

- Harry

"Helps us to establish new business relationships with physicians and surgeons."

More patients, but more importantly, more happy patients and providers is of the utmost importance to us."

- Luis Zuniga, Owner, Fyzical Therapy & Balance Centers, Texas.

"A diagnostic-aid tool that fills the need for objective and dynamic information about the fully loaded, articulating knee."

- Dr Michael Suk,
MD, JD, MPH, MBA, FACS,
Professor and Chair, Geisinger
Musculoskeletal Institute

"KneeKG is very helpful in helping us slow the progression of the knee OA disease, but also treating ACL injuries, anterior knee pain and post operative residual pain."

- Dr Ranjan Sachdev,
MD, FAAOS, MBA, CHC
Member of the International Cartilage
Regeneration and Joint Preservation
Society



"Restoring function via surgery requires much more information than what we currently have. Having more objective information just make sense."

- Dr Michael A. Mont, MD
Board-certified in Orthopedic
Surgery; Fellow, American Academy
of Orthopaedic Surgeons

"The perfect complement to static information derived from X-ray and MRI, providing actionable data on the knee in motion."

- Dr Ian Kurth,
Radiology Associates of Wausau

ENHANCE YOUR PRACTICE NOW

Knee Kinesiography: BETTER INFORMATION. BETTER DECISIONS.**


Report showing results to correct
with ability to compare over time

KNEE KINESIOGRAPHY		
KneeKG Patient ID: 01647 Left knee - Comfortable walking (4.8 km/h)		
Pre and Post Total Knee Replacement (4 months)	Pre-Op 2019-03-25	4 months Post-Op 2019-07-09
Medial compartment		
Varus thrust during loading ^{1,4}	Negative	Negative (Varus 13.4°)
Varus functional lower limb alignment ^{1,4}	Positive (Varus 15.8°)	Positive (Varus 5.4°)
Varus alignment at initial contact ^{1,4}	Positive (Varus 15.8°)	Positive (Varus 6.9°)
Varus alignment during stance ^{1,4}	Positive (Varus 16.3°)	Positive (Varus 7.6°)
Lateral compartment and femoropatellar compartment		
Valgus thrust during loading ^{1,4}	Negative	Negative (Varus 13.4°)
Valgus functional lower limb alignment ^{1,4}	Negative (Varus 13.4°)	Negative (Varus 6.4°)
Valgus alignment at initial contact ^{1,4}	Negative (Varus 15.8°)	Negative (Varus 6.9°)
Valgus alignment during stance ^{1,4}	Negative (Varus 16.3°)	Negative (Varus 7.6°)
Femoropatellar compartment		
Knee flexion at initial contact ^{1,4}	Negative (9.7° of flexion)	Negative (5.1° of flexion)
External tibial rotation at initial contact ^{1,4}	Negative (2.2°)	Negative
General		

Explanations of pain-related
causes and symptoms

APPENDIX: Explanation of the Identified Biomechanical Markers	
Exaggerated extension movement	
General	Possible explanations: <ul style="list-style-type: none"> Enhances patellar instability
Stance	Possible explanations: <ul style="list-style-type: none"> Reveal a lack of calf muscles recruitment which lead to a deficient push-off Possible causes: <ul style="list-style-type: none"> Soleus muscle stiffness or lack of flexibility
Limited flexion movement: Deficient absorption	
General	Possible explanations: <ul style="list-style-type: none"> Limits the efficient impact absorption, thus increasing tibiofemoral and patellofemoral loads Muscular co-contraction requiring an exaggerated and sustained hamstring and quadriceps recruitment Possible causes: <ul style="list-style-type: none"> Proprioception deficit Deficient eccentric quadriceps recruitment

Improved patient management with home
exercise program and objective results



AAOS Guidelines 2021

When integrated in patient education, Knee Kinesiography with the KneeKG system provides statistically significant and clinically superior improvements compared to the results obtained with current standard of care (*Cagnin et al. 2020*)², *as recognized in the practice guidelines published by the American Association of Orthopaedic Surgeons (AAOS), which guidelines are also endorsed by the American Association of Hip and Knee Surgeons (AAHKS)*⁹

About Emovi

We provide actionable intelligence to the orthopedic community. Emovi's KneeKG solution is the first technology to accurately measure the 3D dynamic alignment of the knee. It provides objective data, aids surgeons in their decision-making to correct knee alignment and restore 3D knee function. Our goal is to increase the efficiency and effectiveness of knee joint programs, improve the patient experience, advance population health, enhance the provider experience and contribute to lowering the costs for all stakeholders.¹⁰



*The KneeKG® system is FDA (510k) cleared, Health Canada licensed and CE Marked, to assess the 3D motion of the knee of patients who have impaired movement functions of an orthopaedic cause.

** Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this information does not replace the need for such examination and/or advice in whole or in part. Emovi does not practice medicine. Each physician should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training physicians have received.

- Lustig S, et al. The KneeKG system: a review of the literature. *Knee Surg Sports Traumatol Arthrosc.* 2012 Apr;20(4):633-8.
- Cagnin, A. et al. A multi-arm cluster randomized clinical trial of the use of Knee Kinesiography in the management of osteoarthritis patients in a primary care setting. *Postgrad Med.* 2020 Jan;132(1):91-101.
- D'Souza N, et al. Are biomechanics during gait associated with the structural disease onset and progression of lower limb osteoarthritis? A systematic review and meta-analysis. *Osteoarthritis Cartilage.* 2022 Mar;30(3):381-394.
- Tungtrongjit Y, et al. The effect of preoperative quadriceps exercise on functional outcome after total knee arthroplasty. *J Med Assoc Thai.* 2012 Oct;95 Suppl 10:S58-66.
- 10- Deroche, E., et al. *Arch Orthop Trauma Surg* 142, 1645-1651 (2022)
- Hagemister N, et al. Associations between biomechanical markers from a Knee Kinesiography exam and patient outcome measures: a comprehensive review. *Osteoarthritis & Cartilage.* 2022. V30-S1: S135.

- Planckaert, et al. (2018). Total knee arthroplasty with unexplained pain: new insights from kinematics. *Arch Orthop Trauma Surg.* 138(4): 553-561.
- Larose, et al. (2019). Can total knee arthroplasty restore the correlation between radiographic mechanical axis angle and dynamic coronal plane alignment during gait? *Knee.* 26(3):586-594.
- Wink A.E., et al. (2017). Varus thrust during walking and the risk of incident and worsening medial tibiofemoral MRI lesions: the Multicenter Osteoarthritis Study. *Osteoarthritis and Cartilage.* 25(6): 839-845.
- American Academy of Orthopaedic Surgeons Management of Osteoarthritis of the Knee (NonArthroplasty) Evidence-Based Clinical Practice Guideline. <https://www.aaos.org/pak3cpg> Published 08/31/2021.
- <https://www.milliman.com/en/insight/knee-osteoarthritis-in-a-commercially-insured-population>



Contact us for more information:
info@emovi.ca • 1-866-707-0866
www.emovi.ca